

DAB Enhancer ab675

2 References

概述

产品名称	DAB Enhancer
常规说明	<p>To use as DAB signal enhancer in conjunction with peroxidase based immunostaining systems.</p> <p>i) Add one drop (50µL) of DAB enhancer in 1mL of DAB incubation buffer. ii) Add DAB chromogen and Peroxidase Block (3% Hydrogen Peroxide) to the above buffer. iii) Incubate tissue sections with the incubation medium.</p> <p>DAB is a widely used chromogen for immunoperoxidase staining and immunoblotting. It has been well accepted amongst pathologists because of its superior performance as compared to amino ethylcarbazole (AEC). DAB gives cleaner background as opposed to AEC. Specimens stained in DAB can be dehydrated in ethanol and cleared in xylene and can be mounted for permanent record keeping. Strength of signal generated during the immunostaining is the key to a good staining. This is a stable liquid solution which increases the staining intensity of the tissue specimens stained with DAB several fold and hence increases the efficiency of detection systems.</p>

性能

形式	Liquid
存放说明	Shipped at Room Temperature. Store at Room Temperature.
Reagent说明	<p>DAB is a widely used chromogen for immunoperoxidase staining and immunoblotting. It has been well accepted amongst pathologists because of its superior performance as compared to amino ethylcarbazole (AEC). DAB gives cleaner background as opposed to AEC. Specimens stained in DAB can be dehydrated in ethanol and cleared in xylene and can be mounted for permanent record keeping. Strength of signal generated during the immunostaining is the key to a good staining. This is a stable liquid solution which increases the staining intensity of the tissue specimens stained with DAB several fold and hence increases the efficiency of detection systems.</p>

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

-
- Replacement or refund for products not performing as stated on the datasheet
 - Valid for 12 months from date of delivery
 - Response to your inquiry within 24 hours

 - We provide support in Chinese, English, French, German, Japanese and Spanish
 - Extensive multi-media technical resources to help you
 - We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors