abcam

Product datasheet

Anti-Parainfluenza Virus type 3 antibody ab 28584

1 References

概述

产**品名称** Anti-Parainfluenza Virus type 3抗体

描述 山羊多克隆抗体to Parainfluenza Virus type 3

宿主 Goat

特异性 This antibody detects the parainfluenza virus, type 3. This antibody cross-reacts with type 2, and

has minimal reactivity with type 1.

经测试应用 适用于: ICC/IF

种属反应性 与反应: Human respirovirus 3

免疫原 Tissue, cells or virus corresponding to Human Parainfluenza Virus type 3.

阳性对照 Parainfluenza 3 infected cells.

常规说明

These antibodies assist in research by detecting the infectious disease agent.

The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets

your needs before purchasing.

If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be

found below, along with publications, customer reviews and Q&As

性能

形式 Liquid

存放说明 Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

存储溶液 pH: 7.2

Preservative: 0.1% Sodium azide

Constituent: PBS

纯**度** Protein G purified

Primary antibody说明 These antibodies assist in research by detecting the infectious disease agent.

克隆 多克隆

同种型 IgG

1

The Abpromise guarantee

Abpromise™承诺保证使用ab28584于以下的经测试应用

"应用说明"部分 下显示的仅为推荐的起始稀释度:实际最佳的稀释度/浓度应由使用者检定。

应 用	Ab评论	说明
ICC/IF		Use at an assay dependent concentration.

靶标

相关性

Infectious diseases are the leading cause of death worldwide. AIDS, tuberculosis (TB), malaria, diarrhoeal and respiratory infections account for 78% of deaths caused by infectious disease. As many infectious diseases are controlled, new diseases emerge and old diseases become resistant to current drugs. Many infectious diseases have been associated with an increase risk of carcinoma. Influenza continues to attract researchers as new strains appear by the ability of the influenza gene to mix with different forms of the virus. Recently, research on SARS and West Nile virus has risen due to the increased number of infections.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- · Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit https://www.abcam.cn/abpromise or contact our technical team.

Terms and conditions

• Guarantee only valid for products bought direct from Abcam or one of our authorized distributors