# abcam

### **Product datasheet**

## Human GIPC1 knockout HEK-293T cell line ab266406

1图像

#### 概述

产 <b>品名称</b>	人GIPC1 knockout HEK-293T cell line
Parental Cell Line	HEK293T
Organism	Human
Mutation description	Knockout achieved by using CRISPR/Cas9, Homozygous: Insertion of the selection cassette in exon 5
Passage number	<20
Knockout validation	Sanger Sequencing
Biosafety level	2
常规说明	<b>Recommended control:</b> Human wild-type HEK293T cell line ( <u>ab255449</u> ). Please note a wild- type cell line is not automatically included with a knockout cell line order, if required please add recommended wild-type cell line at no additional cost using the code WILDTYPE-TMTK1.
	<b>Cryopreservation cell medium:</b> Cell Freezing Medium-DMSO Serum free media, contains 8.7% DMSO in MEM supplemented with methyl cellulose.
	Culture medium: DMEM (High Glucose) + 10% FBS
	<b>Initial handling guidelines:</b> Upon arrival, the vial should be stored in liquid nitrogen vapor phase and not at -80°C. Storage at -80°C may result in loss of viability.
	<ol> <li>Thaw the vial in 37°C water bath for approximately 1-2 minutes.</li> <li>Transfer the cell suspension (0.8 mL) to a 15 mL/50 mL conical sterile polypropylene centrifuge tube containing 8.4 mL pre-warmed culture medium, wash vial with an additional 0.8 mL culture medium (total volume 10 mL) to collect remaining cells, and centrifuge at 201 x g (rcf) for 5 minutes at room temperature. 10 mL represents minimum recommended dilution. 20 mL represents maximum recommended dilution.</li> <li>Resuspend the cell pellet in 5 mL pre-warmed culture medium and count using a haemocytometer or alternative cell counting method. Based on cell count, seed cells in an appropriate cell culture flask at a density of 2x10<sup>4</sup> cells/cm<sup>2</sup>. Seeding density is given as a guide only and should be scaled to align with individual lab schedules.</li> <li>Incubate the culture at 37°C incubator with 5% CO<sub>2</sub>. Cultures should be monitored daily.</li> </ol>
	Subculture guidelines:
	All seeding densities should be based on cell counts gained by established methods. A guide seeding density of 2x10 <sup>4</sup> cells/cm <sup>2</sup> is recommended. A partial media change 24 hours prior to subculture may be helpful to encourage growth, if required.

Cells should be passaged when they have achieved 80-90% confluence. This product is subject to limited use licenses from The Broad Institute, ERS Genomics Limited and Sigma-Aldrich Co. LLC, and is developed with patented technology. For full details of the licenses and patents please refer to our **limited use license** and **patent pages**.

We will provide viable cells that proliferate on revival.

#### 性能

Number of cells	1 x 10 <sup>6</sup> cells/vial, 1 mL
Adherent /Suspension	Adherent
Tissue	Kidney
Cell type	epithelial
STR Analysis	Amelogenin X D5S818: 8, 9 D13S317: 12, 14 D7S820: 11 D16S539: 9, 13 vWA: 16, 19 TH01: 7, 9.3 TPOX: 11 CSF1PO: 11, 12
Mycoplasma free	Yes
存 <b>放</b> 说明	Shipped on Dry Ice. Store in liquid nitrogen.
存储溶液	Constituents: 8.7% Dimethylsulfoxide, 2% Cellulose, methyl ether

<b>靶</b> 标	
功能	May be involved in G protein-linked signaling.
组织 <b>特异性</b>	Widely expressed.
序列相似性	Belongs to the GIPC family. Contains 1 PDZ (DHR) domain.
细 <b>胞定位</b>	Cytoplasm. Membrane.

#### 图片



Homozygous: Insertion of the selection cassette in exon5

HEK293T cell line (ab266406)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

#### Our Abpromise to you: Quality guaranteed and expert technical support

- · Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery

- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <u>https://www.abcam.cn/abpromise</u> or contact our technical team.

#### **Terms and conditions**

• Guarantee only valid for products bought direct from Abcam or one of our authorized distributors