

Annexin V/ANXA5-FITC Apoptosis Detection Reagent (500X) ab14082

18 References **1 图像**

概述

产品名称	Annexin V/ANXA5-FITC Apoptosis Detection试剂(500X)
经测试应用	适用于: FM, Flow Cyt, ICC/IF
常规说明	<p>This product is manufactured by BioVision, an Abcam company and was previously called 1001 Annexin V-FITC Reagent.</p> <p>Annexin V/ANXA5-FITC Apoptosis Detection Reagent ab14082 (previously known as Annexin V-FITC Reagent 1001) can be used in a 10 min, one-step staining procedure to detect apoptosis by staining phosphatidylserine molecules which have translocated to the outside of the cell membrane. Analysis is by flow cytometry or fluorescence microscopy.</p> <p>This Annexin V/ANXA5-FITC reagent can be used with Annexin V/ANXA5 binding buffer ab14084 (see protocol section for instructions). The Annexin V/ANXA5-FITC reagent is also available as part of Annexin V/ANXA5-FITC apoptosis kit ab14085.</p> <p>For more apoptosis assays, review the full set of <u>Annexin V/ANXA5 assays</u>, or the <u>apoptosis assay and apoptosis marker guide</u>.</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C. Do Not Freeze.
存储溶液	<p>pH: 7.40</p> <p>Preservative: 0.03% Proclin 300</p> <p>Constituents: 1% BSA, PBS</p>
功能	This protein is an anticoagulant protein that acts as an indirect inhibitor of the thromboplastin-specific complex, which is involved in the blood coagulation cascade.
疾病相关	Pregnancy loss, recurrent, 3
序列相似性	<p>Belongs to the annexin family.</p> <p>Contains 4 annexin repeats.</p>
结构域	<p>The [IL]-x-C-x-x-[DE] motif is a proposed target motif for cysteine S-nitrosylation mediated by the iNOS-S100A8/A9 transnitrosylase complex.</p> <p>A pair of annexin repeats may form one binding site for calcium and phospholipid.</p>
翻译后修饰	S-nitrosylation is induced by interferon-gamma and oxidatively-modified low-density lipoprotein

(LDL(ox)) possibly implicating the iNOS-S100A8/9 transnitrosylase complex.

应用


The Abpromise guarantee

Abpromise™ 承诺保证使用 ab14082 于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
FM		Use at an assay dependent concentration.
Flow Cyt		Use at an assay dependent dilution.
ICC/IF		Use 1µl for 10 ⁵ cells.

图片



Murine platelets were stimulated with CRP (5 µg/ml) in the presence of annexin V-fluorescein isothiocyanate for 10 minutes. The annexin V-positive population is indicated (AnV +ve).

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors