

Anti-Enterovirus 71 antibody [10F0] ab36367

18 References

概述	
产品名称	Anti-Enterovirus 71抗体[10F0]
描述	小鼠单克隆抗体[10F0] to Enterovirus 71
宿主	Mouse
经测试应用	适用于: ICC/IF, Neutralising
种属反应性	与反应: Enterovirus A71
免疫原	Tissue, cells or virus corresponding to Enterovirus 71. Native Enterovirus 71
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能	
形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
存储溶液	Constituent: Ascites
纯度	Tissue culture supernatant
克隆	单克隆
克隆编号	10F0
同种型	IgG1

应用	
The Abpromise guarantee	<u>Abpromise™</u> 承诺保证使用ab36367于以下的经测试应用
“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。	

应用	Ab评论	说明
ICC/IF		1/500 - 1/2000. Fix with cold acetone for 10 minutes. Dilute antibody in PBS (pH 7.5 to pH 8.0) with 1% BSA.
Neutralising		1/14.

靶标

相关性 Enterovirus 71 was first isolated in 1969. Since then it has been associated with sporadic cases of outbreaks in various parts of the world, including the U.S., Brazil, Europe, Australia, Malaysia, and Taiwan. Enterovirus 71 infection may present itself much like other enterovirus infections, with symptoms ranging from asymptomatic to the presentation of rashes, diarrhea, colds, meningitis, encephalitis, hand to foot to mouth disease, myocarditis, or some combination of these syndromes. On the other hand, Enterovirus 71 may cause many unusual syndromes such as acute flaccid paralysis (similar to paralytic poliomyelitis), bulbar and brain encephalitis, Guillain Barre syndrome, and rapidly fatal pulmonary edema and hemorrhage.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors