

Product datasheet

Anti-6X His tag® antibody (Agarose) ab1231

概述

产品名称	Anti-6X His tag®抗体(Agarose)
描述	山羊多克隆抗体to 6X His tag® (Agarose)
宿主	Goat
偶联物	Agarose
经测试应用	适用于: IP 不适用于: ChIP
免疫原	Full length native HHHHHH (6-His)(purified): conjugated to KLH. Run BLAST with Run BLAST with
常规说明	HIS-TAG® is a trademark of EMD Biosciences, Inc. Affinity purified antibodies were coupled to agarose beads using a cyanogen bromide method.

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	Preservative: 0.1% Sodium Azide Constituents: 0.1M Sodium chloride, 0.01M PBS, pH 7.2
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab1231** in the following tested applications.
The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
----	------	----

应用	Ab评论	说明
IP		Use at an assay dependent concentration. Use 20-40 µl of gel slurry per 0.1-1 mg of protein lysate or extract.
应用说明		Is unsuitable for ChIP.
靶标		
相关性		The H-H-H-H-H-H motif is used as a tag on many recombinant proteins to facilitate purification. The antibody recognizes the His-tag fused to the amino- or carboxy- termini of targeted proteins in transfected or transformed cells.
细胞定位		Depends upon the localization of the parent protein tagged with hexahistidine.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors