


Product datasheet

Anti-ABLIM2 antibody ab100926

2 图像

概述

产品名称	Anti-ABLIM2抗体
描述	兔多克隆抗体to ABLIM2
宿主	Rabbit
经测试应用	适用于: WB, IHC-P
种属反应性	与反应: Human 预测可用于: Mouse, Rat 
免疫原	Recombinant fragment, corresponding to amino acids 74-309 of Human ABLIM2 (BC122567).
阳性对照	Human Fetal Heart, Human Fetal Brain, Human Fetal Muscle, Human Fetal Skeletal Muscle

性能

形式	Lyophilised:Reconstitute with 200ul distilled sterile water. Please note that if you receive this product in liquid form it has already been reconstituted as described and no further reconstitution is necessary.
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
存储溶液	Preservative: 0.02% Sodium Azide Constituents: PBS, pH 7.2
纯度	Protein A purified
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab100926** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		1/500 - 1/1000. Predicted molecular weight: 68 kDa.
IHC-P		1/100 - 1/500.

靶标

功能	May act as scaffold protein. May stimulate ABRA activity and ABRA-dependent SRF transcriptional activity.
组织特异性	Highly expressed in skeletal muscle.
序列相似性	Contains 1 HP (headpiece) domain. Contains 4 LIM zinc-binding domains.
细胞定位	Cytoplasm. In skeletal muscle, sarcomeric or cosarcomeric localization.

图片




All lanes : Anti-ABLIM2 antibody (ab100926)
at 1/500 dilution

Lane 1 : Human fetal heart lysate

Lane 2 : Human fetal brain lysate

Lane 3 : Human fetal muscle lysate

Predicted band size: 68 kDa



Cytoplasmic staining of ABLIM2 in a Formalin/PFA-fixed paraffin-embedded section of Human Fetal Skeletal Muscle using ab100926 at a dilution of 1/100.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-ABLIM2 antibody (ab100926)

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish

- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors