

## Product datasheet

# Anti-ALPK2 antibody ab111909

### 1 图像

#### 概述

产品名称	Anti-ALPK2抗体
描述	兔多克隆抗体to ALPK2
宿主	Rabbit
经测试应用	适用于: WB, IHC-P
种属反应性	与反应: Human
免疫原	A peptide from C terminal residues of Human ALPK2 protein (NM_052947.3).
阳性对照	Human fetal liver lysate.

#### 性能

形式	Lyophilised:Reconstitute in 200ul Sterile Distilled Water.
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.
存储溶液	pH: 7.20 Preservative: 0.02% Sodium azide Constituents: 1% BSA, 98.98% PBS
纯度	Immunogen affinity purified
纯化说明	ab111909 is purified by a peptide affinity column.
克隆	多克隆
同种型	IgG

#### 应用

Our [Abpromise guarantee](#) covers the use of **ab111909** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		1/500 - 1/1000. Detects a band of approximately 237 kDa (predicted molecular weight: 237 kDa).
IHC-P		1/100 - 1/500.

## 靶标

功能	Kinases that recognize phosphorylation sites in which the surrounding peptides have an alpha-helical conformation.
序列相似性	Belongs to the protein kinase superfamily. Alpha-type protein kinase family. ALPK subfamily. Contains 1 alpha-type protein kinase domain. Contains 2 Ig-like (immunoglobulin-like) domains.

## 图片



Anti-ALPK2 antibody (ab111909) at 1/500 dilution + Human fetal liver lysate

**Predicted band size:** 237 kDa

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors