


Product datasheet

# Anti-Apolipoprotein H antibody ab115215

1 图像

概述

产品名称	Anti-Apolipoprotein H抗体
描述	山羊多克隆抗体to Apolipoprotein H
宿主	Goat
特异性	Recognizes human Apolipoprotein H. A single precipitin arc is observed with normal plasma. No arc was seen with Apolipoprotein H deficient plasma.
经测试应用	适用于: WB, IHC-P
种属反应性	与反应: Human 预测可用于: Mouse, Rat, Rabbit, Dog, Pig 
免疫原	Human Apolipoprotein H purified from plasma.
阳性对照	Human Liver tissue

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.
存储溶液	pH: 7.40 Constituents: 0.24% HEPES, 0.88% Sodium chloride, 50% Glycerol
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab115215** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


应用	Ab评论	说明
WB		Use at an assay dependent dilution. Predicted molecular weight: 38 kDa.

应用	Ab评论	说明
IHC-P		Use a concentration of 2.5 µg/ml. Perform heat mediated antigen retrieval with citrate buffer pH 6 before commencing with IHC staining protocol.

## 靶标

功能	Binds to various kinds of negatively charged substances such as heparin, phospholipids, and dextran sulfate. May prevent activation of the intrinsic blood coagulation cascade by binding to phospholipids on the surface of damaged cells.
组织特异性	Expressed by the liver and secreted in plasma.
序列相似性	Contains 4 Sushi (CCP/SCR) domains.
翻译后修饰	N- and O-glycosylated. PubMed:6587378 also reports glycosylation on 'Asn-188' for their allele.
细胞定位	Secreted.

## 图片



Human liver (formalin fixed, paraffin embedded) stained with ab115215 at 2.5 ug/ml followed by biotinylated anti Goat IgG secondary antibody, alkaline phosphatase streptavidin and chromogen.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-Apolipoprotein H antibody (ab115215)

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors