

## Product datasheet

# Anti-Atrophin 1 antibody ab104727

### 1 图像

#### 概述

<b>产品名称</b>	Anti-Atrophin 1抗体
<b>描述</b>	兔多克隆抗体to Atrophin 1
<b>经测试应用</b>	<b>适用于:</b> WB
<b>种属反应性</b>	<b>与反应:</b> Human <b>预测可用于:</b> Mouse, Rat, Rabbit, Horse, Cow, Cat, Dog, Pig, Chimpanzee
<b>免疫原</b>	Synthetic peptide corresponding to a region within N terminal amino acids 38-87 (VSTSSSDGKA EKSRQTAKKA RVEEASTPKV NKQGRSEEIS ESESEETNAP) of Human Atrophin 1 (NP_001931). <a href="#">Run BLAST with ExPASy</a> <a href="#">Run BLAST with NCBI</a>
<b>阳性对照</b>	Placenta lysate

#### 性能

<b>形式</b>	Liquid
<b>存放说明</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.
<b>存储溶液</b>	Preservative: 0.09% Sodium azide Constituents: PBS, 2% Sucrose
<b>纯度</b>	Immunogen affinity purified
<b>克隆</b>	多克隆
<b>同种型</b>	IgG

#### 应用

Our [Abpromise guarantee](#) covers the use of **ab104727** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		Use a concentration of 1 µg/ml. Predicted molecular weight: 125 kDa. Good results were obtained when blocked with 5% non-fat dry milk in 0.05% PBS-T.

## 靶标

相关性	Atrophin 1 (ATN1) is essential for the development of multiple tissues. Defects in ATN1 are the cause of dentatorubral-pallidoluysian atrophy (DRPLA).
细胞定位	Cytoplasm. Nuclear

## 图片



Anti-Atrophin 1 antibody (ab104727) at 1  $\mu$ g/ml + placenta lysate at 10  $\mu$ g

**Predicted band size : 125 kDa**

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors