

Product datasheet

Anti-BAG2 antibody [EPR3568] ab80596

重组 RabMAb

2 图像

概述

产品名称	Anti-BAG2抗体[EPR3568]
描述	兔单克隆抗体[EPR3568] to BAG2
经测试应用	适用于: WB, ICC/IF 不适用于: Flow Cyt, IHC-P or IP
种属反应性	与反应: Mouse, Rat, Human
免疫原	Synthetic peptide from the C terminal of Human BAG2
阳性对照	WB: Jurkat, HeLa and A549 cell lysate; ICC/IF: HeLa cells.
常规说明	This product is a recombinant rabbit monoclonal antibody. Our RabMAb [®] technology is a patented hybridoma-based technology for making rabbit monoclonal antibodies. For details on our patents, please refer to RabMAb[®] patents

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at -20°C. Stable for 12 months at -20°C.
存储溶液	PBS 49%, Sodium azide 0.01%, Glycerol 50%, BSA 0.05%
纯度	Tissue culture supernatant
克隆	单克隆
克隆编号	EPR3568
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab80596** in the following tested applications.
The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		1/1000 - 1/2000. Detects a band of approximately 24 kDa (predicted molecular weight: 24 kDa).
ICC/IF		1/100 - 1/250.


应用说明 Is unsuitable for Flow Cyt, IHC-P or IP.

靶标

功能 Inhibits the chaperone activity of HSP70/HSC70 by promoting substrate release.


序列相似性 Contains 1 BAG domain.

图片



Immunofluorescence analysis of HeLa cells with 1/100 ab80596.

Immunocytochemistry/ Immunofluorescence - BAG2 antibody [EPR3568] (ab80596)



All lanes : Anti-BAG2 antibody [EPR3568] (ab80596) at 1/2000 dilution

Lane 1 : Jurkat cell lysate

Lane 2 : HeLa cell lysate

Lane 3 : A549 cell lysate

Lysates/proteins at 10 µg per lane.

Secondary

goat anti-rabbit HRP at 1/2000 dilution

Predicted band size : 24 kDa

Observed band size : 24 kDa

Western blot - BAG2 antibody [EPR3568] (ab80596)

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors