


Product datasheet

Anti-C13orf31 antibody ab87684

1 图像

概述

产品名称	Anti-C13orf31抗体
描述	兔多克隆抗体to C13orf31
经测试应用	适用于: WB
种属反应性	与反应: Human 预测可用于: Mouse, Rat, Rabbit, Horse, Guinea pig, Cow, Cat, Dog, Pig 
免疫原	Synthetic peptide corresponding to a region within the internal amino acids 180-229 (IITSSLIPDI FHHGFTTRTG GISYIPTLSS FNLFSSSKRR DPKVVVQENL) of Human C13orf31 (NP_694950). Run BLAST with ExPASy Run BLAST with NCBI
阳性对照	Fetal lung lysate.

性能


形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.
存储溶液	Preservative: None Constituents: 2% Sucrose, PBS
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab87684** in the following tested applications.
The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab 评论	说明
WB		Use a concentration of 1 µg/ml. Detects a band of approximately 48 kDa (predicted molecular weight: 48 kDa). Good results were obtained when blocked with 5% non-fat dry milk in 0.05% PBS-T.

Anti-C13orf31 antibody 图像



Anti-C13orf31 antibody (ab87684) at 1 $\mu\text{g/ml}$
(in 5% skim milk / PBS buffer) + human fetal
lung lysate at 10 μg

Secondary

HRP conjugated anti-Rabbit IgG at 1/50000
dilution

Predicted band size : 48 kDa

Observed band size : 48 kDa

Additional bands at : ~30 kDa. We are
unsure as to the identity of these extra bands.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors