


Product datasheet

Anti-C17orf76 antibody ab121998

2 图像

概述

产品名称	Anti-C17orf76抗体
描述	兔多克隆抗体to C17orf76
经测试应用	适用于: ICC/IF, IHC-P
种属反应性	与反应: Human 预测可用于: Mouse 
免疫原	Recombinant fragment corresponding to amino acids 56-125 (GMVQELLRMV RQGRREEAGT LLQHLRQDLG MESTSLDDVL YRYASFRNLV DPITHDLIIS LARYIHCPKP) of Human C17orf76 (Q7TSF4). Run BLAST with ExPASy Run BLAST with NCBI
阳性对照	Human lateral ventricle tissue.

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.
存储溶液	pH: 7.20 Preservative: 0.02% Sodium azide Constituents: 59% PBS, 40% Glycerol
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab121998** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明


应用	Ab评论	说明
ICC/IF		Use a concentration of 1 - 4 µg/ml. Recommend PFA Fixation and Triton X-100 treatment
IHC-P		1/20 - 1/50. Perform heat mediated antigen retrieval with citrate buffer pH 6 before commencing with IHC staining protocol. pH6

靶标

序列相似性


Contains 2 LRR (leucine-rich) repeats.

图片



Immunofluorescent staining of Human cell line A-431 shows positivity in nucleus but not nucleoli and vesicles. Recommended concentration of ab121998 1-4 µg/ml. Cells treated with PFA/Triton X-100.

Immunocytochemistry/ Immunofluorescence - Anti-C17orf76 antibody (ab121998)



ab121998, at a 1/25 dilution, staining C17orf76 in paraffin-embedded Human lateral ventricle tissue by Immunohistochemistry.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-C17orf76 antibody (ab121998)

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors