


Product datasheet

Anti-C1orf177 antibody ab87303

1 图像

概述

产品名称	Anti-C1orf177抗体
描述	兔多克隆抗体to C1orf177
经测试应用	适用于: WB
种属反应性	与反应: Human 预测可用于: Rabbit, Horse, Cow 
免疫原	Synthetic peptide, corresponding to a region within internal sequence amino acids 252-301 9SMQKKKPRELMNFKSFVEELNSHHNKKHGVSFKLPRNPKTPTERIWAN L) of Human C1orf177 (NP_689820)
阳性对照	721_B cell lysate.

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.
存储溶液	Preservative: None Constituents: 2% Sucrose, PBS
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab87303** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		Use a concentration of 1 µg/ml. Predicted molecular weight: 47 kDa. Good results were obtained when blocked with 5% non-fat dry milk in 0.05% PBS-T.

## 靶标


---

## 相关性

The specific function of C1orf177 is not yet known.

## 图片

---



Anti-C1orf177 antibody (ab87303) at 1 µg/ml  
+ 721\_B cell lysate at 10 µg

### Secondary

anti-Rabbit IgG HRP at 1/50000 dilution

**Predicted band size :** 47 kDa

**Observed band size :** 47 kDa

**Additional bands at :** 80 kDa. We are unsure as to the identity of these extra bands.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors