

Product datasheet

Anti-C6orf62 antibody ab107371

3 图像

概述

产品名称	Anti-C6orf62抗体
描述	兔多克隆抗体to C6orf62
经测试应用	适用于: WB, IHC-P
种属反应性	与反应: Mouse, Hamster, Human
免疫原	Synthetic peptide conjugated to KLH, between amino acids 118-147 from the central region of Human C6orf62 (NP_112201.1).
阳性对照	CHO cell line lysate, Mouse kidney tissue lysate, Human brain tissue

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at 4°C (up to 6 months). Store at -20°C long term.
存储溶液	Preservative: 0.09% Sodium Azide Constituents: PBS
纯度	Immunogen affinity purified
纯化说明	ab107371 is purified through a protein A column, followed by peptide affinity purification.
克隆	多克隆
同种型	IgG

应用


Our [Abpromise guarantee](#) covers the use of **ab107371** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		1/100 - 1/500. Predicted molecular weight: 27 kDa.
IHC-P		1/50 - 1/100. Perform heat mediated antigen retrieval with citrate buffer pH 6 before commencing with IHC staining protocol.


靶标

图片




Anti-C6orf62 antibody (ab107371) at 1/100 dilution + CHO cell line lysate at 35 µg

Predicted band size : 27 kDa



Anti-C6orf62 antibody (ab107371) at 1/100 dilution + Mouse kidney tissue lysate at 35 µg

Predicted band size : 27 kDa



Immunohistochemistry analysis in formalin fixed and paraffin embedded Human brain tissue followed by peroxidase conjugation of the secondary antibody and DAB staining. ab107371 at a dilution of 1/50.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery

- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors