

Product datasheet

Anti-CCDC75 antibody ab122451

3 图像

概述

产品名称	Anti-CCDC75抗体
描述	兔多克隆抗体to CCDC75
经测试应用	适用于: WB, IHC-P, ICC/IF
种属反应性	与反应: Mouse, Rat, Human
免疫原	antigen sequence: KSGIGHEASL KRKAEEKLES YRKKIHMKNQ AEEKAAEQFR MRLKNKQDEM KLEGDLRRSQ RACQQLDVQK NIQVPREAWY WL, corresponding to amino acids 107-188 of Human CCDC75 Run BLAST with ExPASy Run BLAST with NCBI
阳性对照	Human testis tissue; RT4, U251 MG, Human Liver and Tonsil lysates

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.
存储溶液	pH: 7.20 Preservative: 0.02% Sodium azide Constituents: 59% PBS, 40% Glycerol
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab122451** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		1/250 - 1/500.


应用	Ab评论	说明
IHC-P		1/200 - 1/500. Perform heat mediated antigen retrieval with citrate buffer pH 6 before commencing with IHC staining protocol.
ICC/IF		Use a concentration of 1 - 4 µg/ml. Recommend PFA Fixation and Triton X-100 treatment

靶标

序列相似性


Contains 1 G-patch domain.

图片




Immunofluorescent staining of Human cell line U-2 OS shows positivity in nucleus.
Recommended concentration of ab122451 1-4 µg/ml. Cells treated with PFA/Triton X-100.

Immunocytochemistry/ Immunofluorescence - Anti-CCDC75 antibody (ab122451)



ab122451, at 1/200, staining CCDC75 in paraffin embedded Human testis tissue by immunohistochemistry.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-CCDC75 antibody (ab122451)



Western blot - Anti-CCDC75 antibody (ab122451)

All lanes : Anti-CCDC75 antibody
(ab122451) at 1/250 dilution

Lane 1 : RT4 lysate

Lane 2 : U251 MG lysate

Lane 3 : Human plasma lysate

Lane 4 : Liver lysate

Lane 5 : Tonsil lysate

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors