

Product datasheet

CCRF-CEM nuclear extract lysate ab14642

概述

产品名称	CCRF-CEM核extract裂解物
常规说明	We recommend aliquoting the extracts into single use fractions and then storing.

性能

形式	Liquid
存放说明	Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles.
存储溶液	Preservative: None Constituents: 20% Glycerol, 20mM HEPES, 100mM Potassium Chloride, 1mM Magnesium Chloride, 0.5mM PMSF, 0.5mM DTT. pH 7.9
背景	The direct isolation and propagation of a continuous culture of human lymphoblasts (CCRF-CEM) (also called CEM) from the peripheral blood of a patient with acute lymphoblastic leukemia has permitted extensive sequential cytogenetic analysis of these cells during continuous propagation in culture.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors