

Product datasheet

Anti-CD11a antibody [B-B15] (Phycoerythrin) ab27351

1 References

概述

产品名称	Anti-CD11a抗体[B-B15] (Phycoerythrin)
描述	小鼠单克隆抗体[B-B15] to CD11a (Phycoerythrin)
偶联物	Phycoerythrin. Ex: 488nm, Em: 575nm
经测试应用	适用于: Flow Cyt
种属反应性	与反应: Human
免疫原	Tissue/ cell preparation (Human), HL60 cell line.

性能

形式	Lyophilised
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	Preservative: 0.1% Sodium Azide Constituents: 1% BSA, PBS
纯度	Ion Exchange Chromatography
克隆	单克隆
克隆编号	B-B15
骨髓瘤	x63-Ag8.653
同种型	IgG1

应用

Our [Abpromise guarantee](#) covers the use of **ab27351** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
----	------	----

Flow Cyt

应用说明	Flow Cyt: Use 10 µl to label 10 <sup>6</sup> cells or 100 µl of whole blood.  Not yet tested in other applications.
------	---------------------------------------------------------------------------------------------------------------------------

Optimal dilutions/concentrations should be determined by the end user.

## 靶标

---

<b>功能</b>	Integrin alpha-L/beta-2 is a receptor for ICAM1, ICAM2, ICAM3 and ICAM4. It is involved in a variety of immune phenomena including leukocyte-endothelial cell interaction, cytotoxic T-cell mediated killing, and antibody dependent killing by granulocytes and monocytes.
<b>组织特异性</b>	Leukocytes.
<b>序列相似性</b>	Belongs to the integrin alpha chain family. Contains 7 FG-GAP repeats. Contains 1 VWFA domain.
<b>结构域</b>	The integrin I-domain (insert) is a VWFA domain. Integrins with I-domains do not undergo protease cleavage.
<b>细胞定位</b>	Membrane.

---

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors