

## Product datasheet

# Anti-CD134 / OX40L receptor antibody [OX-86] (FITC) ab33998

### 概述

产品名称	Anti-CD134 / OX40L receptor抗体[OX-86] (FITC)
描述	大鼠单克隆抗体[OX-86] to CD134 / OX40L receptor (FITC)
宿主	Rat
偶联物	FITC. Ex: 493nm, Em: 528nm
经测试应用	适用于: Flow Cyt
种属反应性	与反应: Mouse
免疫原	Tissue, cells or virus corresponding to Mouse CD134/ OX40L receptor. CHO cells transfected with murine CD134/ OX40L receptor.

### 性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	Preservative: 0.09% Sodium Azide Constituents: 1.0% BSA, PBS, pH 7.4
纯度	Protein G purified
克隆	单克隆
克隆编号	OX-86
骨髓瘤	NS1
同种型	IgG1

### 应用

Our [Abpromise guarantee](#) covers the use of **ab33998** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
Flow Cyt		Use 10µl for 10 <sup>6</sup> cells. <a href="#">ab18404</a> - Rat monoclonal IgG1, is suitable for use as an isotype control with this antibody.

## 靶标

---

功能	Receptor for TNFSF4/OX40L/GP34.
序列相似性	Contains 4 TNFR-Cys repeats.
细胞定位	Membrane.

---

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors