

## Product datasheet

# Anti-CD134 / OX40L receptor antibody [OX-86] (FITC) ab33998

### 概述

|       |   |
|-------|---|
| 产品名称  | Anti-CD134 / OX40L receptor抗体[OX-86] (FITC)   |
| 描述    | 大鼠单克隆抗体[OX-86] to CD134 / OX40L receptor (FITC)   |
| 宿主    | Rat   |
| 偶联物   | FITC. Ex: 493nm, Em: 528nm  |
| 经测试应用 | 适用于: Flow Cyt   |
| 种属反应性 | 与反应: Mouse  |
| 免疫原   | Tissue, cells or virus corresponding to Mouse CD134/ OX40L receptor. CHO cells transfected with murine CD134/ OX40L receptor. |

### 性能

|      |   |
|------|---|
| 形式   | Liquid  |
| 存放说明 | Shipped at 4°C. Store at +4°C.  |
| 存储溶液 | Preservative: 0.09% Sodium Azide<br>Constituents: 1.0% BSA, PBS, pH 7.4 |
| 纯度   | Protein G purified  |
| 克隆   | 单克隆   |
| 克隆编号 | OX-86   |
| 骨髓瘤  | NS1   |
| 同种型  | IgG1  |

### 应用

Our [Abpromise guarantee](#) covers the use of **ab33998** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

| 应用       | Ab评论 | 说明  |
|----------|------|---|
| Flow Cyt |      | Use 10µl for 10 <sup>6</sup> cells.<br><a href="#">ab18404</a> - Rat monoclonal IgG1, is suitable for use as an isotype control with this antibody. |

## 靶标

---

|       |                                 |
|-------|---------------------------------|
| 功能    | Receptor for TNFSF4/OX40L/GP34. |
| 序列相似性 | Contains 4 TNFR-Cys repeats.    |
| 细胞定位  | Membrane.                       |

---

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors