

Product datasheet

Anti-CD28 antibody [B-T3] (FITC) ab27310

1 References 1 图像

概述

产品名称	Anti-CD28抗体[B-T3] (FITC)
描述	小鼠单克隆抗体[B-T3] to CD28 (FITC)
偶联物	FITC. Ex: 493nm, Em: 528nm
特异性	Recognises the B7/BB1 receptor (B7/BB1 R)
经测试应用	适用于: Flow Cyt
种属反应性	与反应: Human
免疫原	U266 cell line

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	Preservative: 0.1% Sodium Azide Constituents: 5% BSA, PBS
纯度	Ion Exchange Chromatography
克隆	单克隆
克隆编号	B-T3
骨髓瘤	x63-Ag8.653
同种型	IgG2a

应用

Our [Abpromise guarantee](#) covers the use of **ab27310** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
----	------	----

Flow Cyt

应用说明 Flow Cyt: Use 10 µl to label 10⁶ cells or 100 µl of whole blood.


Not yet tested in other applications.

Optimal dilutions/concentrations should be determined by the end user.

靶标

功能	Involved in T-cell activation, the induction of cell proliferation and cytokine production and promotion of T-cell survival.
组织特异性	Expressed in T-cells and plasma cells, but not in less mature B-cells.
序列相似性	Contains 1 Ig-like V-type (immunoglobulin-like) domain.
细胞定位	Membrane.

图片



A typical staining pattern with ab27310 of lymphocytes

Flow Cytometry / FACS - CD28 antibody [B-T3]
(FITC) (ab27310)

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors