

Product datasheet

Anti-CD5 antibody [B-B8] ab27371

1 References 1 图像

概述

产品名称	Anti-CD5抗体[B-B8]
描述	小鼠单克隆抗体[B-B8] to CD5
宿主	Mouse
特异性	Recognises CD5.
经测试应用	适用于: Flow Cyt
种属反应性	与反应: Human
免疫原	PHA activated T cells.

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	Preservative: 0.1% Sodium Azide Constituents: 1% BSA, PBS
纯度	Ion Exchange Chromatography
克隆	单克隆
克隆编号	B-B8
骨髓瘤	x63-Ag8.653
同种型	IgG1

应用

Our [Abpromise guarantee](#) covers the use of **ab27371** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
----	------	----

Flow Cyt

应用说明 Flow Cyt: Use 10µl for 10⁶ cells or 100µl of whole blood.


Not yet tested in other applications.

Optimal dilutions/concentrations should be determined by the end user.

靶标

功能	May act as a receptor in regulating T-cell proliferation. CD5 interacts with CD72/LYB-2.
序列相似性	Contains 3 SRCR domains.
细胞定位	Membrane.

图片



Ab27371, at a dilution of 10 μ l for 10⁶ cells or 100 μ l of whole blood, staining CD5 by Flow Cytometry.

Flow Cytometry - Anti-CD5 antibody [B-B8]
(ab27371)

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors