

Product datasheet

Anti-CDC42EP3 antibody - N-terminal ab93594

1 References 1 图像

概述

| | |
|-------|--|
| 产品名称 | Anti-CDC42EP3抗体- N-terminal |
| 描述 | 兔多克隆抗体to CDC42EP3 - N-terminal |
| 宿主 | Rabbit |
| 经测试应用 | 适用于: WB, ELISA |
| 种属反应性 | 与反应: Human |
| 免疫原 | Synthetic peptide selected from the N terminal region of Human CDC42EP3 conjugated to KLH (NP_006440.2). |
| 阳性对照 | 293 cell line lysates. |

性能

| | |
|------|--|
| 形式 | Liquid |
| 存放说明 | Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles. |
| 存储溶液 | Preservative: 0.09% Sodium Azide Constituents: PBS |
| 纯度 | Immunogen affinity purified |
| 纯化说明 | ab93594 is purified through a protein A column, followed by peptide affinity purification. |
| 克隆 | 多克隆 |
| 同种型 | IgG |

应用

Our [Abpromise guarantee](#) covers the use of **ab93594** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

| 应用 | Ab评论 | 说明 |
|-------|------|----|
| WB | | |
| ELISA | | |

应用说明

ELISA: 1/1000.

WB: 1/100 - 1/500. Predicted molecular weight: 28 kDa.

Not yet tested in other applications.

Optimal dilutions/concentrations should be determined by the end user.

靶标


相关性

CDC42EP3 is probably involved in the organization of the actin cytoskeleton. May act downstream of CDC42 to induce actin filament assembly leading to cell shape changes. It induces pseudopodia formation in fibroblasts.

细胞定位

Endomembrane system; Peripheral membrane protein. Cytoplasm; cytoskeleton.

图片



Anti-CDC42EP3 antibody - N-terminal
(ab93594) at 1/100 dilution + 293 cell line
lysates at 35 µg

Predicted band size: 28 kDa

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors