


Product datasheet

Anti-COQ5 antibody ab122936

1 图像

概述

产品名称	Anti-COQ5抗体
描述	兔多克隆抗体to COQ5
经测试应用	适用于: WB
种属反应性	与反应: Rat 预测可用于: Mouse, Rabbit, Horse, Cow, Cat, Dog, Human, Pig 
免疫原	Synthetic peptide corresponding to a region within N-terminal amino acids 41-90 (RSLSQEKRAAETHFGFETVSEKEKGGKVYQVFQSVARKYDLMNDMMSLG I) of Rat COQ5 (NP_001034111). Run BLAST with Run BLAST with
阳性对照	Rat Lung lysate

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
存储溶液	Constituents: 97% PBS, 2% Sucrose
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab122936** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		Use a concentration of 1 µg/ml. Predicted molecular weight: 37 kDa. Good results were obtained when blocked with 5% non-fat dry milk in 0.05% PBS-T.

靶标

功能	Methyltransferase required for the conversion of 2-polyprenyl-6-methoxy-1,4-benzoquinol (DDMQH2) to 2-polyprenyl-3-methyl-6-methoxy-1,4-benzoquinol (DMQH2).
通路	Cofactor biosynthesis; ubiquinone biosynthesis.
序列相似性	Belongs to the methyltransferase superfamily. UbiE family.
细胞定位	Mitochondrion.

图片



Anti-COQ5 antibody (ab122936) at 1 μ g/ml +

Rat Lung lysate at 10 μ g

Predicted band size : 37 kDa

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors