

## Product datasheet

# Anti-CRIM1 antibody [MM0202-16G1] ab89293

### 概述

产品名称	Anti-CRIM1抗体[MM0202-16G1]
描述	小鼠单克隆抗体[MM0202-16G1] to CRIM1
宿主	Mouse
经测试应用	适用于: WB
种属反应性	与反应: Human
免疫原	Recombinant Human CRIM1 extracellular domain.
常规说明	Reconstituted antibody should be aliquoted and can be stored frozen at < -20°C for at least six months without detectable loss of activity. Avoid repeated freeze-thaw cycles.

### 性能

形式	Lyophilised:Reconstitute the antibody with 200 µl sterile PBS to give a final concentration of 0.5mg/ml.
存放说明	Shipped at 4°C. Store at -80°C.
存储溶液	Preservative: None Constituents: PBS
纯度	Protein G purified
纯化说明	The IgG fraction of culture supernatant was purified by Protein G affinity chromatography.
克隆	单克隆
克隆编号	MM0202-16G1
同种型	IgG1

### 应用

Our [Abpromise guarantee](#) covers the use of **ab89293** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		1/500 - 1/1000. Predicted molecular weight: 114 kDa.

## 靶标

---

<b>功能</b>	May play a role in CNS development by interacting with growth factors implicated in motor neuron differentiation and survival. May play a role in capillary formation and maintenance during angiogenesis. Modulates BMP activity by affecting its processing and delivery to the cell surface.
<b>组织特异性</b>	Expressed in pancreas, kidney, skeletal muscle, lung, placenta, brain, heart, spleen, liver and small intestine. Expressed in blood vessels (at protein level).
<b>序列相似性</b>	Contains 4 antistasin-like domains. Contains 1 IGFBP N-terminal domain. Contains 6 VWFC domains.
<b>翻译后修饰</b>	N-glycosylated.
<b>细胞定位</b>	Secreted and Cell membrane.

---

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors