

Product datasheet

Anti-Complement factor B antibody [M13/12] ab106139

2 References

概述

产品名称	Anti-Complement factor B抗体[M13/12]
描述	小鼠单克隆抗体[M13/12] to Complement factor B
宿主	Mouse
特异性	ab106139 recognizes complement factor B/Bb
经测试应用	适用于: WB, IP
种属反应性	与反应: Human
免疫原	Complement factor B/Bb

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
存储溶液	Preservative: 0.02% Sodium Azide Constituents: 0.1% BSA, PBS
纯度	Protein G purified
克隆	单克隆
克隆编号	M13/12
同种型	IgG1

应用

Our [Abpromise guarantee](#) covers the use of **ab106139** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		1/50. Predicted molecular weight: 86 kDa.
IP		Use at an assay dependent dilution.

靶标

相关性	Factor B which is part of the alternate pathway of the complement system is cleaved by factor D into 2 fragments: Ba and Bb. Bb, a serine protease, then combines with complement factor 3b to generate the C3 or C5 convertase. It has also been implicated in proliferation and differentiation of preactivated B-lymphocytes, rapid spreading of peripheral blood monocytes, stimulation of lymphocyte blastogenesis and lysis of erythrocytes. Ba inhibits the proliferation of preactivated B-lymphocytes.
细胞定位	Plasma membrane and Secreted

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors