

Product datasheet

Anti-Complement factor B antibody [M13/12] ab106139

2 References

概述

| | |
|-------|--|
| 产品名称 | Anti-Complement factor B抗体[M13/12] |
| 描述 | 小鼠单克隆抗体[M13/12] to Complement factor B |
| 宿主 | Mouse |
| 特异性 | ab106139 recognizes complement factor B/Bb |
| 经测试应用 | 适用于: WB, IP |
| 种属反应性 | 与反应: Human |
| 免疫原 | Complement factor B/Bb |

性能

| | |
|------|---|
| 形式 | Liquid |
| 存放说明 | Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles. |
| 存储溶液 | Preservative: 0.02% Sodium Azide Constituents: 0.1% BSA, PBS |
| 纯度 | Protein G purified |
| 克隆 | 单克隆 |
| 克隆编号 | M13/12 |
| 同种型 | IgG1 |

应用

Our [Abpromise guarantee](#) covers the use of **ab106139** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

| 应用 | Ab评论 | 说明 |
|----|------|---|
| WB | | 1/50. Predicted molecular weight: 86 kDa. |
| IP | | Use at an assay dependent dilution. |

靶标

| | |
|-------------|---|
| 相关性 | Factor B which is part of the alternate pathway of the complement system is cleaved by factor D into 2 fragments: Ba and Bb. Bb, a serine protease, then combines with complement factor 3b to generate the C3 or C5 convertase. It has also been implicated in proliferation and differentiation of preactivated B-lymphocytes, rapid spreading of peripheral blood monocytes, stimulation of lymphocyte blastogenesis and lysis of erythrocytes. Ba inhibits the proliferation of preactivated B-lymphocytes. |
| 细胞定位 | Plasma membrane and Secreted |

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors