

Product datasheet

Anti-Cytomegalovirus pp65 antibody [CH12] ab53489

2 References

概述

| | |
|-------|--|
| 产品名称 | Anti-Cytomegalovirus pp65抗体[CH12] |
| 描述 | 小鼠单克隆抗体[CH12] to Cytomegalovirus pp65 |
| 宿主 | Mouse |
| 特异性 | ab53489 is reactive with pp65 (UL83) of Cytomegalovirus. |
| 经测试应用 | 适用于: ICC/IF, WB |
| 种属反应性 | Reacts with Cytomegalovirus. |
| 免疫原 | HCMV AD169 IC Extract. |

性能

| | |
|------|---|
| 形式 | Liquid |
| 存放说明 | Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles. |
| 存储溶液 | Preservative: None Constituents: PBS, pH 7.4 |
| 纯度 | Protein G purified |
| 克隆 | 单克隆 |
| 克隆编号 | CH12 |
| 同种型 | IgG |

应用

Our [Abpromise guarantee](#) covers the use of **ab53489** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

| 应用 | Ab评论 | 说明 |
|--------|------|----------------------------------|
| ICC/IF | | Use a concentration of 10 µg/ml. |
| WB | | Use a concentration of 10 µg/ml. |

靶标

相关性

Cytomegalovirus is a member of the herpes virus group, which includes herpes simplex virus types 1 and 2, varicella zoster virus (which causes chicken pox), and Epstein Barr virus (which causes infectious mononucleosis). These viruses share a characteristic ability to remain dormant within the body over a long period. CMV viral genes are co-ordinately expressed in groups at various times after infection. Early viral proteins are expressed in the nucleus of infected cells within 3 to 24 hours of infection prior to the commencement of viral DNA replication. This is followed by expression of the early intermediate genes, which encode enzymes required for viral DNA replication. After 48 to 72 hours, a number of late viral antigens may be demonstrated in the nuclei and cytoplasm of infected cells. pp65 is a 65kD phosphorylated glycoprotein and is the most abundant of the late antigens.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors