


Product datasheet

Anti-DDX26B antibody ab95083

1 References 1 图像

概述

产品名称	Anti-DDX26B抗体
描述	兔多克隆抗体to DDX26B
经测试应用	适用于: WB
种属反应性	与反应: Human 预测可用于: Mouse, Rat, Rabbit, Horse, Chicken, Guinea pig, Cow, Cat, Dog, Zebrafish 
免疫原	Synthetic peptide corresponding to a region within the N terminal amino acids 180-229 (STEPEQLGSV PTDESAITQM CEVTGGRSYC VRTQRMLNQC LESLVQKVQS) of Human DDX26B (NP_872346). Run BLAST with ExPASy Run BLAST with NCBI
阳性对照	THP1 cell lysate

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
存储溶液	Preservative: None Constituents: 2% Sucrose, PBS
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab95083** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		Use a concentration of 1 µg/ml. Predicted molecular weight: 97 kDa. Good results were obtained when blocked with 5% non-fat dry milk in 0.05% PBS-T.

相关性

The function of this protein remains unknown. There are 2 isoforms produced by alternative splicing.

Anti-DDX26B antibody 图像



Anti-DDX26B antibody (ab95083) at 1 μ g/ml
(in 5% skim milk / PBS buffer) + THP1 cell
lysate at 10 μ g

Secondary

HRP conjugated anti-Rabbit IgG at 1/50000
dilution

Predicted band size : 97 kDa

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors