

Product datasheet

Anti-DGKK antibody ab111042

2 图像

概述

产品名称	Anti-DGKK抗体
描述	兔多克隆抗体to DGKK
宿主	Rabbit
经测试应用	适用于: IHC-P, ICC/IF
种属反应性	与反应: Human
免疫原	Synthetic peptide derived from the C-terminal of Human DGKK.
阳性对照	COS7 cells; Human brain tissue.

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at -20°C. Stable for 12 months at -20°C.
存储溶液	pH: 7.40 Preservative: 0.02% Sodium azide Constituents: 49.1% PBS, 0.88% Sodium chloride, 50% Glycerol Note: without Mg ²⁺ and Ca ²⁺
纯度	Immunogen affinity purified
纯化说明	ab111042 was affinity-purified from Rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab111042** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明


应用	Ab评论	说明
IHC-P		1/50 - 1/100. Perform heat mediated antigen retrieval with citrate buffer pH 6 before commencing with IHC staining protocol.
ICC/IF		1/100 - 1/500.

靶标

相关性 DGKK (Diacylglycerol kinase kappa) phosphorylates diacylglycerol (DAG) to generate phosphatidic acid (PA).

细胞定位 Cytoplasm. Membrane; Peripheral membrane protein

图片



ab111042 at 1/100 dilution staining DGKK in COS7 cells by Immunofluorescence. The image on the right is treated with the synthesized peptide.



ab111042 at 1/50 dilution staining DGKK in paraffin-embedded Human brain tissue by Immunohistochemistry. The image on the right is treated with the synthesized peptide.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Promise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish

- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors