

## Product datasheet

# Anti-DMRT1 antibody ab105954

### 1 图像

#### 概述

产品名称	Anti-DMRT1抗体
描述	兔多克隆抗体to DMRT1
宿主	Rabbit
经测试应用	适用于: WB
种属反应性	与反应: Mouse 预测可用于: Rat
免疫原	Synthetic peptide corresponding to a region within the C terminal amino acids 321-370 (VFSPSSQDS GLVSLSSSSP MSNESSKGVL ECESASSEPS SYAVNQVLEE) of Mouse DMRT1 (NP_056641). <a href="#">Run BLAST with ExPASy</a> <a href="#">Run BLAST with NCBI</a>

#### 性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.
存储溶液	Preservative: None Constituents: 2% Sucrose, PBS
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

#### 应用

Our [Abpromise guarantee](#) covers the use of **ab105954** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		Use a concentration of 1 µg/ml. Predicted molecular weight: 39 kDa. Good results were obtained when blocked with 5% non-fat dry milk in 0.05% PBS-T.

## 靶标

功能	May be required for testis development.
组织特异性	Testis specific.
序列相似性	Belongs to the DMRT family. Contains 1 DM DNA-binding domain.
发展阶段	Expressed in the genital ridges of six-week-old male embryos. Becomes confined to the forming seminiferous tubules (possibly Sertoli cells) at seven weeks. Not detected in female embryos.
细胞定位	Nucleus.

## 图片



Anti-DMRT1 antibody (ab105954) (0.2-1 ug) +  
NIH/3T3 cell lysate

**Predicted band size: 39 kDa**

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors