

Product datasheet

Anti-DnaJ antibody ab80160

1 References

概述

产品名称	Anti-DnaJ抗体
描述	兔多克隆抗体to DnaJ
宿主	Rabbit
经测试应用	适用于: WB
种属反应性	与反应: Chlamydomonas reinhardtii
免疫原	Recombinant DnaJ of Chlamydomonas reinhardtii

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at -20°C. Stable for 12 months at -20°C.
存储溶液	Preservative: None Constituents: Whole serum
纯度	Whole antiserum
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab80160** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		1/5000. Predicted molecular weight: 41 kDa.

靶标

相关性 The *E. coli* heat-shock protein DnaJ has been implicated in protein folding and protein complex dissociation. DnaJ, Heat shock protein, functions in association with DnaK(Hsp70) molecular

chaperone to facilitate protein folding. p70 chaperone. DnaJ plays a key role in the chaperone reaction by stimulating the ATPase activity and activating the substrate binding of Hsp70. DnaJ consists of four domains that are N-terminal 76 amino acid J-domain, G/F domain, zinc-binding cysteine rich CR-domain, C-terminal CTD-domain and they are conserved to various degrees among the homologues.

细胞定位

Cytoplasmic

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors