

## Product datasheet

# Donkey Anti-Human IgG H&L (HRP) preadsorbed ab102429

### 概述

|               |   |
|---------------|---|
| <b>产品名称</b>   | 驴抗人IgG H&L (HRP)预吸附二抗   |
| <b>宿主</b>     | Donkey  |
| <b>靶标种属</b>   | Human   |
| <b>特异性</b>    | By immunoelectrophoresis and ELISA this antibody reacts specifically with human IgG and with light chains common to other human immunoglobulins. No antibody was detected against non-immunoglobulin serum proteins. Less than 1% cross reactivity to bovine, chicken, goat, horse, mouse, rabbit, rat and sheep IgG was detected. This antibody may cross react with IgG from other species. |
| <b>经测试应用</b>  | <b>适用于:</b> IHC-P, ICC, ELISA, WB   |
| <b>最小交叉反应</b> | Chicken, Cow, Goat, Horse, Mouse, Rabbit, Rat, Sheep <a href="#">more details</a>   |
| <b>偶联物</b>    | HRP   |

### 性能

|             |   |
|-------------|---|
| <b>形式</b>   | Liquid  |
| <b>存放说明</b> | Shipped at 4°C. Store at +4°C.  |
| <b>存储溶液</b> | Preservative: 0.1% Proclin<br>Constituents: 0.2% BSA, PBS   |
| <b>纯度</b>   | Immunogen affinity purified   |
| <b>纯化说明</b> | Antiserum was cross adsorbed using bovine, chicken, goat, horse, mouse, rabbit, rat and sheep immunosorbents to remove cross reactive antibodies. This antibody was isolated by affinity chromatography using antigen coupled to agarose beads and conjugated to HRP. |
| <b>克隆</b>   | 多克隆   |
| <b>同种型</b>  | IgG   |

### 应用

Our [Abpromise guarantee](#) covers the use of **ab102429** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

| 应用    | Ab评论 | 说明   |
|-------|------|--|
| IHC-P |      | 1/200 - 1/5000.  |
| ICC   |      | 1/200 - 1/5000.  |
| ELISA |      | 1/10000 - 1/100000. Primary ELISA.   |
| WB    |      | 1/5000 - 1/50000. Colorimetric: 1/5000 - 1/30000, Chemiluminescent: 1/10000 - 1/50000. |

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

### Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors