

Product datasheet

Donkey Anti-Human IgG H&L (HRP) preadsorbed ab102429

概述

产品名称	驴Anti-Human IgG H&L (HRP) preadsorbed
靶标种属	Human
特异性	By immunoelectrophoresis and ELISA this antibody reacts specifically with human IgG and with light chains common to other human immunoglobulins. No antibody was detected against non-immunoglobulin serum proteins. Less than 1% cross reactivity to bovine, chicken, goat, horse, mouse, rabbit, rat and sheep IgG was detected. This antibody may cross react with IgG from other species.
经测试应用	适用于: IHC-P, ICC, ELISA, WB
最小交叉反应	Chicken, Cow, Goat, Horse, Mouse, Rabbit, Rat, Sheep more details
偶联物	HRP

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	Preservative: 0.1% Proclin Constituents: 0.2% BSA, PBS
纯度	Immunogen affinity purified
纯化说明	Antiserum was cross adsorbed using bovine, chicken, goat, horse, mouse, rabbit, rat and sheep immunosorbents to remove cross reactive antibodies. This antibody was isolated by affinity chromatography using antigen coupled to agarose beads and conjugated to HRP.
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab102429** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
IHC-P		1/200 - 1/5000.

应用	Ab评论	说明
ICC		1/200 - 1/5000.
ELISA		1/10000 - 1/100000. Primary ELISA.
WB		1/5000 - 1/50000. Colorimetric: 1/5000 - 1/30000, Chemiluminescent: 1/10000 - 1/50000.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors