

## Product datasheet

# Donkey Anti-Rabbit IgG H&L (Alkaline Phosphatase) preadsorbed ab97082

### 1 References

#### 概述

<b>产品名称</b>	驴抗兔IgG H&L (Alkaline Phosphatase)预吸附二抗
<b>描述</b>	驴多克隆抗体二抗to兔IgG - H&L (AP),预吸附二抗
<b>宿主</b>	Donkey
<b>靶标种属</b>	Rabbit
<b>特异性</b>	By immunoelectrophoresis and ELISA this antibody reacts specifically with Rabbit IgG and with light chains common to other Rabbit immunoglobulins. No antibody was detected against non immunoglobulin serum proteins. Less than 0.1% cross reactivity to bovine, chicken, goat, horse, human, mouse, pig and rat was detected.
<b>经测试应用</b>	<b>适用于:</b> ICC, WB, IHC-P, ELISA
<b>最小交叉反应</b>	Chicken, Cow, Goat, Human, Mouse, Pig, Rat <a href="#">more details</a>
<b>免疫原</b>	Other Immunogen Type corresponding to Rabbit IgG.
<b>偶联物</b>	Alkaline Phosphatase

#### 性能

<b>形式</b>	Liquid
<b>存放说明</b>	Shipped at 4°C. Store at +4°C.
<b>存储溶液</b>	Preservative: 0.09% Sodium Azide Constituents: 0.2% BSA, 0.1M Sodium chloride, 50mM HEPES, 1mM Magnesium chloride, 0.1mM Zinc chloride, pH 7.1
<b>纯度</b>	Immunogen affinity purified
<b>纯化说明</b>	Antiserum was cross adsorbed using bovine, chicken, goat, human, mouse, pig and rat immunosorbents to remove cross reactive antibodies. This antibody was isolated by affinity chromatography using antigen coupled to agarose beads and conjugated to Alkaline Phosphatase (AP).
<b>克隆</b>	多克隆
<b>同种型</b>	IgG

#### 应用

Our [Abpromise guarantee](#) covers the use of **ab97082** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
ICC		Use at an assay dependent dilution.
WB		1/2500 - 1/25000. (Chemiluminescent)
IHC-P		1/100 - 1/1000.
ELISA		1/2500 - 1/25000. (Primary ELISA)

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

### Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors