

## Product datasheet

# Donkey Anti-Sheep IgG H&L (Alkaline Phosphatase) ab6901

### 概述

产品名称	驴Anti-Sheep IgG H&L (Alkaline Phosphatase)
靶标种属	Sheep
经测试应用	适用于: ICC/IF, IHC-Fr, IHC-P, WB, ELISA, Immunomicroscopy, Dot blot
免疫原	Full length native Sheep IgG (purified).
偶联物	Alkaline Phosphatase

### 性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	Preservative: 0.01% Sodium Azide Constituents: 50% Glycerol, 10mg/ml BSA, 0.15M Sodium chloride, 0.05M Tris HCl, 1mM Magnesium chloride, 0.1mM Zinc chloride, pH 8
纯度	IgG fraction
纯化说明	This product was prepared from monospecific antiserum by immunoaffinity chromatography using Sheep IgG coupled to agarose beads followed by solid phase adsorption(s) to remove any unwanted reactivities
共轭说明	Alkaline Phosphatase (Calf Intestine) (Molecular Weight 140,000 daltons)
克隆	多克隆

### 应用

Our [Abpromise guarantee](#) covers the use of **ab6901** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
ICC/IF		
IHC-Fr		
IHC-P		

应用	Ab评论	说明
WB		
ELISA		
Immunomicroscopy		
Dot blot		

#### 应用说明

Dot: Use at an assay dependant dilution.

ELISA: 1/3000 - 1/14000.

IHC-P: Use at an assay dependant dilution.

IHC-Fr: Use at an assay dependant dilution.

IM: Use at an assay dependant dilution.

WB: Use at an assay dependant dilution.

Not tested in other applications.

Optimal dilutions/concentrations should be determined by the end user.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

#### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

#### Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors