

Product datasheet

Donkey Anti-Sheep IgG H&L (DyLight® 594) ab96941

1 图像

概述

产品名称	驴抗羊IgG H&L (DyLight® 594)
描述	驴多克隆抗体二抗to羊IgG - H&L (DyLight® 594)
宿主	Donkey
靶标种属	Sheep
特异性	By immunoelectrophoresis and ELISA this antibody reacts specifically with Sheep IgG and with light chains common to other Sheep immunoglobulins. No antibody was detected against non immunoglobulin serum proteins.
经测试应用	适用于: IHC-P, ICC/IF, Flow Cyt
偶联物	DyLight® 594. Ex: 593nm, Em: 618nm

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	Preservative: 0.09% Sodium Azide Constituents: 0.2% BSA, PBS
纯度	Immunogen affinity purified
纯化说明	This antibody was isolated by affinity chromatography using antigen coupled to agarose beads and conjugated to DyLight® 594.
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab96941** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


应用	Ab评论	说明
IHC-P		1/50 - 1/500.
ICC/IF		1/50 - 1/500.

应用	Ab评论	说明
----	------	----

Flow Cyt

1/50 - 1/200.

图片



Emission spectra of DyLight® fluorochromes available in our catalog.

Line colors represent the approximate visible colors of the wavelength maxima.

DyLight®-Donkey polyclonal Secondary Antibody to Sheep IgG - H&L (DyLight® 594)(ab96941)

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors