


Product datasheet

Anti-E2 230K antibody ab118832

1 图像

概述

产品名称	Anti-E2 230K抗体
描述	山羊多克隆抗体to E2 230K
经测试应用	适用于: ELISA, IHC-P
种属反应性	与反应: Human 预测可用于: Mouse, Rat, Cow, Dog, Pig 
免疫原	Synthetic peptide: ALLEAGMPECTEDK , corresponding to C terminal amino acids 1279-1292 of Human E2 230K (NP_071349.2). Run BLAST with Run BLAST with
阳性对照	Human skeletal muscle tissue.

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.
存储溶液	pH: 7.30 Preservative: 0.02% Sodium azide Constituents: 99% Tris buffered saline, 0.5% BSA
纯度	Protein G purified
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab118832** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


应用	Ab评论	说明
ELISA		1/32000.

应用	Ab评论	说明
IHC-P		Use a concentration of 3.75 µg/ml.

靶标

功能	Catalyzes the covalent attachment of ubiquitin to other proteins.
组织特异性	Predominantly expressed in skeletal muscle and heart.
通路	Protein modification; protein ubiquitination.
序列相似性	Belongs to the ubiquitin-conjugating enzyme family.

图片



ab118832, at 3.75 µg/ml, staining E2 230K in formalin fixed, paraffin embedded Human skeletal muscle tissue by Immunohistochemistry.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-E2 230K antibody (ab118832)

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors