


Product datasheet

Anti-ECHDC1 antibody ab124449

2 图像

概述

产品名称	Anti-ECHDC1抗体
描述	兔多克隆抗体to ECHDC1
宿主	Rabbit
经测试应用	适用于: WB, IHC-P
种属反应性	与反应: Mouse, Human 预测可用于: Rat 
免疫原	Recombinant fragment, corresponding to amino acids 74-269 of Human ECHDC1 (BC003549).
阳性对照	Human fetal liver and kidney lysates, Mouse liver lysate, Human fetal skin tissue

性能

形式	Lyophilised:Reconstitute in 200ul sterile H2O.
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.
存储溶液	pH: 7.20 Preservative: 0.02% Sodium azide Constituents: 98% PBS, 1% BSA
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab124449** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


应用	Ab评论	说明
WB		1/1000 - 1/2000. Predicted molecular weight: 34 kDa.
IHC-P		1/100 - 1/500.

靶标

序列相似性

Belongs to the enoyl-CoA hydratase/isomerase family.

图片




All lanes : Anti-ECHDC1 antibody
(ab124449) at 1/1000 dilution

Lane 1 : Human fetal liver lysate

Lane 2 : Mouse liver lysate

Lane 3 : Human fetal kidney lysate

Predicted band size: 34 kDa



ab124449 at 1/100 dilution staining ECHDC1
in Human fetal skin tissue by
Immunohistochemistry Formalin-fixed,
Paraffin-embedded tissue.

Immunohistochemistry (Formalin/PFA-fixed paraffin-
embedded sections) - Anti-ECHDC1 antibody
(ab124449)

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors