

Product datasheet

Anti-EFCAB1 antibody ab121758

2 图像

概述

产品名称	Anti-EFCAB1抗体
描述	兔多克隆抗体to EFCAB1
宿主	Rabbit
经测试应用	适用于: IHC-P, WB
种属反应性	与反应: Human
免疫原	antigen sequence, corresponding to amino acids 107-206 (YCFEVDLNG DGFISKEEMF HMLKNSLLKQ PSEEDPDEGI) of Human EFCAB1. <a href="#">Run BLAST with ExPASy</a> <a href="#">Run BLAST with NCBI</a>
阳性对照	EFCAB1-transfected HEK293T lysate, Human duodenum tissue.

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
存储溶液	pH: 7.20 Preservative: 0.02% Sodium azide Constituents: 59% PBS, 40% Glycerol
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab121758** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
IHC-P		1/10 - 1/20. Perform heat mediated antigen retrieval with citrate buffer pH 6 before commencing with IHC staining protocol.


应用	Ab评论	说明
----	------	----

WB 1/250 - 1/500.

## 靶标

**序列相似性** Contains 3 EF-hand domains.

## 图片




**All lanes** : Anti-EFCAB1 antibody (ab121758) at 1/250 dilution

**Lane 1** : vector only transfected HEK293T lysate

**Lane 2** : EFCAB1-transfected HEK293T lysate (Co-expressed with a C-terminal myc-DDK tag (~3.1 kDa))

Developed using the ECL technique.

Western blot - Anti-EFCAB1 antibody (ab121758)



[ab1211758](#) at 1/15 dilution staining EFCAB1 in Paraffin-embedded Human duodenum tissue by Immunohistochemistry.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-EFCAB1 antibody (ab121758)

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery

- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors