

Product datasheet

Anti-EHBP1L1 antibody ab122557

2 图像

概述

|       |  |
|-------|--|
| 产品名称  | Anti-EHBP1L1抗体   |
| 描述    | 兔多克隆抗体to EHBP1L1   |
| 经测试应用 | 适用于: ICC/IF, IHC-P   |
| 种属反应性 | 与反应: Human   |
| 免疫原   | antigen sequence: IEVGLLGVLG IETGAAEGAI LGTQEIASRD SGVPGLEADT TGIQVKEVGG SEVPEIATGT AETEILGTQE IASRSSGVPG LESEVAGAQE TEVGGSGISG, corresponding to amino acids of Human EHBP1L1.<br><a href="#">Run BLAST with ExPASy</a> <a href="#">Run BLAST with NCBI</a> |
| 阳性对照  | Human stomach tissue.  |

性能

|      |   |
|------|---|
| 形式   | Liquid  |
| 存放说明 | Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles. |
| 存储溶液 | pH: 7.20<br>Preservative: 0.02% Sodium azide<br>Constituents: 59% PBS, 40% Glycerol   |
| 纯度   | Immunogen affinity purified   |
| 克隆   | 多克隆   |
| 同种型  | IgG   |

应用

Our [Abpromise guarantee](#) covers the use of **ab122557** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

| 应用     | Ab评论 | 说明   |
|--------|------|--|
| ICC/IF |      | Use a concentration of 1 - 4 µg/ml.<br>Recommend PFA Fixation and Triton X-100 treatment |


| 应用    | Ab评论 | 说明  |
|-------|------|---|
| IHC-P |      | 1/2500 - 1/5000. Perform heat mediated antigen retrieval with citrate buffer pH 6 before commencing with IHC staining protocol. |

## 靶标

### 序列相似性


Contains 1 CH (calponin-homology) domain.

## 图片



Immunofluorescent staining of Human cell line A-431 shows positivity in nucleoli, plasma membrane and cytoplasm. Recommended concentration of ab122557 1-4  $\mu\text{g/ml}$ . Cells treated with PFA/Triton X-100.

Immunocytochemistry/ Immunofluorescence - Anti-EHBP1L1 antibody (ab122557)



ab122557, at 1/2500 dilution, staining EHBP1L1 in paraffin-embedded Human stomach tissue by Immunohistochemistry.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-EHBP1L1 antibody (ab122557)

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

## **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors