

## Product datasheet

### EHD4 peptide (532-541) (Carboxyterminal end) ab45789

#### 概述

---

产品名称 EHD4肽(532-541) (Carboxyterminal end)

#### 描述

---

性质 Synthetic

#### 氨基酸序列

种属 Human

序列 C-SHRKSLPKAD

氨基酸 532 to 541

#### 技术指标

---

Our [Abpromise guarantee](#) covers the use of **ab45789** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用 Blocking

形式 Liquid

#### 制备和贮存

---

稳定性和存储 Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

#### 常规信息

---

相关性 EHD4 is involved in the control of trafficking at the early endosome and regulates exit of cargo toward both the recycling compartment and the late endocytic pathway.

细胞定位 Early endosome membrane; Peripheral membrane protein. Recycling endosome membrane; Peripheral membrane protein.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

## **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors