

## Product datasheet

# Anti-EME1 antibody ab88878

★☆☆☆☆ 2 Abreviews 3 图像

### 概述

产品名称	Anti-EME1抗体
描述	小鼠多克隆抗体to EME1
宿主	Mouse
经测试应用	适用于: WB, ICC/IF
种属反应性	与反应: Human
免疫原	Full length protein, corresponding to amino acids 1 - 583 of Human EME1 (AAH16470.1)
阳性对照	HeLa or Jurkat cell lysate. HeLa cells.

### 性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
存储溶液	Preservative: None Constituents: 1X PBS, pH 7.2
纯度	Protein A purified
克隆	多克隆
同种型	IgG

### 应用

Our [Abpromise guarantee](#) covers the use of **ab88878** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB	★☆☆☆☆	Use a concentration of 1 µg/ml. Predicted molecular weight: 63 kDa.
ICC/IF		Use a concentration of 10 µg/ml.

## 靶标

### 功能

Interacts with MUS81 to form a DNA structure-specific endonuclease with substrate preference for branched DNA structures with a 5'-end at the branch nick. Typical substrates include 3'-flap structures, replication forks and nicked Holliday junctions. May be required in mitosis for the processing of stalled or collapsed replication forks.


### 序列相似性

Belongs to the EME1/MMS4 family.

### 细胞定位


Nucleus > nucleolus. Recruited to regions of DNA damage in S-phase cells.

## 图片



Anti-EME1 antibody (ab88878) at 1 µg/ml +  
HeLa cell lysate at 50 µg

**Predicted band size:** 63 kDa



Anti-EME1 antibody (ab88878) at 1 µg/ml +  
Jurkat cell lysate  
at 50 µg

**Predicted band size:** 63 kDa



ab88878 at 10 µg/ml staining EME1 in HeLa  
cells by immunofluorescence.

### **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors