

## Product datasheet

# Anti-FAM26F antibody ab98871

1 References 1 图像

### 概述

产品名称	Anti-FAM26F抗体
描述	兔多克隆抗体to FAM26F
经测试应用	适用于: WB
种属反应性	与反应: Human 预测可用于: Rabbit
免疫原	Synthetic peptide corresponding to a region within internal amino acids 143-192 (NRSCAAELPL VPCNQAKASD VQDLLKDLKA QSQVLGWILI AVVILLIF) of Human FAM26F (NP_001010919). <a href="#">Run BLAST with ExPASy</a> <a href="#">Run BLAST with NCBI</a>
阳性对照	Human fetal muscle lysate

### 性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
存储溶液	Preservative: None Constituents: 2% Sucrose, PBS
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

### 应用

Our [Abpromise guarantee](#) covers the use of **ab98871** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		Use a concentration of 1 µg/ml. Predicted molecular weight: 35 kDa. Good results were obtained when blocked with 5% non-fat dry milk in 0.05% PBS-T.

## 靶标

功能	Pore-forming subunit of a voltage-gated ion channel.
序列相似性	Belongs to the FAM26 family.
细胞定位	Membrane.

## 图片



Anti-FAM26F antibody (ab98871) at 1  $\mu$ g/ml +  
Human fetal muscle lysate at 10  $\mu$ g

**Predicted band size : 35 kDa**

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors