

Product datasheet

Anti-FBXO36 antibody ab81317

1 图像

概述

产品名称	Anti-FBXO36抗体
描述	兔多克隆抗体to FBXO36
经测试应用	适用于: WB, ELISA
种属反应性	与反应: Human 预测可用于: Mouse, Rat, Rabbit, Horse, Cow, Cat, Pig 
免疫原	Synthetic peptide corresponding to a region within N terminal amino acids 37-86 (WKISLRSEYRSTKPGGEAKETHEDFLENSHLQGQTALIFGARILDYVINLC) of Human FBXO36 (NP_777559). Run BLAST with Run BLAST with
阳性对照	Human fetal liver lysate.

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.
存储溶液	Preservative: None Constituents: 2% Sucrose, PBS
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab81317** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab 评论	说明

应用	Ab 评论	说明
----	-------	----

WB Use a concentration of 1 µg/ml. Predicted molecular weight: 22 kDa. Can be blocked with [FBXO36 peptide \(ab128001\)](#). Good results were obtained when blocked with 5% non-fat dry milk in 0.05% PBS-T.


ELISA Use at an assay dependent concentration.
Titre using peptide based assay: 1/312500.

靶标

功能 Substrate-recognition component of the SCF (SKP1-CUL1-F-box protein)-type E3 ubiquitin ligase complex.

序列相似性 Contains 1 F-box domain.

图片



Anti-FBXO36 antibody (ab81317) at 1 µg/ml
(in 5% skim milk / PBS buffer) + Human fetal liver lysate at 10 µg

Secondary

HRP conjugated anti-Rabbit IgG at 1/50000 dilution

Predicted band size : 22 kDa

Observed band size : 27 kDa

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors