


Product datasheet

Anti-FBXO36 antibody ab81317

1 图像

概述

<b>产品名称</b>	Anti-FBXO36抗体
<b>描述</b>	兔多克隆抗体to FBXO36
<b>经测试应用</b>	<b>适用于:</b> WB, ELISA
<b>种属反应性</b>	<b>与反应:</b> Human <b>预测可用于:</b> Mouse, Rat, Rabbit, Horse, Cow, Cat, Pig 
<b>免疫原</b>	Synthetic peptide corresponding to a region within N terminal amino acids 37-86 (WKISLRSEYRSTKPGGEAKETHEDFLENSHLQGQTALIFGARILDYVINLC) of Human FBXO36 (NP_777559). <a href="#">Run BLAST with</a> <a href="#">Run BLAST with</a>
<b>阳性对照</b>	Human fetal liver lysate.

性能

<b>形式</b>	Liquid
<b>存放说明</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.
<b>存储溶液</b>	Preservative: None Constituents: 2% Sucrose, PBS
<b>纯度</b>	Immunogen affinity purified
<b>克隆</b>	多克隆
<b>同种型</b>	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab81317** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab 评论	说明

应用	Ab 评论	说明
----	-------	----

WB Use a concentration of 1 µg/ml. Predicted molecular weight: 22 kDa. Can be blocked with [FBXO36 peptide \(ab128001\)](#). Good results were obtained when blocked with 5% non-fat dry milk in 0.05% PBS-T.


ELISA Use at an assay dependent concentration.  
Titre using peptide based assay: 1/312500.

## 靶标

**功能** Substrate-recognition component of the SCF (SKP1-CUL1-F-box protein)-type E3 ubiquitin ligase complex.

**序列相似性** Contains 1 F-box domain.

## 图片



Anti-FBXO36 antibody (ab81317) at 1 µg/ml  
(in 5% skim milk / PBS buffer) + Human fetal liver lysate at 10 µg

### Secondary

HRP conjugated anti-Rabbit IgG at 1/50000 dilution

**Predicted band size** : 22 kDa

**Observed band size** : 27 kDa

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors