

Product datasheet

Anti-FPGT antibody ab97914

2 图像

概述

产品名称	Anti-FPGT抗体
描述	兔多克隆抗体to FPGT
宿主	Rabbit
经测试应用	适用于: WB, IHC-P
种属反应性	与反应: Human
免疫原	Recombinant protein fragment, corresponding to a region within amino acids 146-423 of Human FPGT (NP_003829).
阳性对照	293T, A431, H1299, HeLaS3 or Raji cell lines.

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
存储溶液	Preservative: 0.01% Thimerosal (merthiolate) Constituents: 10% Glycerol, 0.1M Tris, 0.1M Glycine, pH 7.0
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab97914** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		1/500 - 1/3000. Predicted molecular weight: 67 kDa.
IHC-P		1/100.

## 靶标

功能	Catalyzes the formation of GDP-L-fucose from GTP and L-fucose-1-phosphate. Functions as a salvage pathway to reutilize L-fucose arising from the turnover of glycoproteins and glycolipids.
组织特异性	Expressed in many tissues.
细胞定位	Cytoplasm.


## 图片



Anti-FPGT antibody (ab97914) at 1/500 dilution + Raji whole cell lysate at 30  $\mu$ g

**Predicted band size:** 67 kDa

7.5% SDS PAGE



ab97914, at 1/100 dilution, staining FPGT in MDA-MB-468 xenograft tissue by Immunohistochemistry. Formalin-fixed, paraffin-embedded section.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-FPGT antibody (ab97914)

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish

- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors