


Product datasheet

Anti-FRG1 antibody ab111616

1 图像

概述

产品名称	Anti-FRG1抗体
描述	兔多克隆抗体to FRG1
经测试应用	适用于: WB
种属反应性	与反应: Mouse 预测可用于: Rat, Rabbit, Goat, Horse, Chicken, Guinea pig, Cow, Cat, Dog, Zebrafish 
免疫原	Synthetic peptide corresponding to a region within internal amino acids 118-167 (GINS DGLVVG RSDAIGPREQ WEPVFQDGKM ALLASNSCFI RCNEAGDIEA) of Mouse FRG1 (NP_038550). <a href="#">Run BLAST with ExPASy</a> <a href="#">Run BLAST with NCBI</a>
阳性对照	Mouse brain lysate.

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.
存储溶液	Constituents: 98% PBS, 2% Sucrose
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab111616** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		Use a concentration of 0.2 - 1 µg/ml. Predicted molecular weight: 29 kDa. Good results were obtained when blocked with 5% non-fat dry milk in 0.05% PBS-T.

## 靶标

<b>功能</b>	May have a role in processing of pre-rRNA or in the assembly of rRNA into ribosomal subunits. May be involved in pre-mRNA splicing.
<b>组织特异性</b>	Expressed in adult muscle, lymphocytes, fetal brain, muscle, and placenta.
<b>序列相似性</b>	Belongs to the FRG1 family.
<b>细胞定位</b>	Nucleus > Cajal body. Nucleus speckle. Nucleus > nucleolus. First expressed in Cajal bodies and nuclear speckles. As speckle expression diminishes, expression is seen in dense structures (dense fibrillar component and the granular component) of the nucleolus.

## 图片



Anti-FRG1 antibody (ab111616) at 1 µg/ml +  
Mouse brain lysate at 10 µg

**Predicted band size : 29 kDa**

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors