

Product datasheet

Anti-FRRS1 antibody ab121538

2 图像

概述

产品名称	Anti-FRRS1抗体
描述	兔多克隆抗体to FRRS1
宿主	Rabbit
经测试应用	适用于: IHC-P, WB
种属反应性	与反应: Human
免疫原	antigen sequence, corresponding to internal sequence amino acids 187-336 of Human FRRS1.
阳性对照	Human epididymis tissue. RT4 cell and U-251 MG cells lysate; Human Plasma and Liver tissue lysate.

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
存储溶液	pH: 7.20 Preservative: 0.02% Sodium azide Constituents: 59% PBS, 40% Glycerol
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab121538** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


应用	Ab评论	说明
IHC-P		1/50 - 1/200. Perform heat mediated antigen retrieval with citrate buffer pH 6 before commencing with IHC staining protocol.

应用	Ab评论	说明
WB		1/250 - 1/500.

## 靶标


<b>功能</b>	Ferric-chelate reductases reduce Fe(3+) to Fe(2+) before its transport from the endosome to the cytoplasm.
<b>序列相似性</b>	Belongs to the FRRS1 family. Contains 1 cytochrome b561 domain. Contains 1 DOMON domain. Contains 1 reelin domain.
<b>细胞定位</b>	Membrane.

## 图片



ab121538 at 1/100 dilution staining FRRS1 in paraffin-embedded Human epididymis tissue by Immunohistochemistry.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-FRRS1 antibody (ab121538)



**All lanes** : Anti-FRRS1 antibody (ab121538) at 1/250 dilution

**Lane 1** : RT4 cell lysate

**Lane 2** : U-251 MG cell lysate

**Lane 3** : Human Plasma lysate

**Lane 4** : Human Liver tissue lysate

**Lane 5** : Human Tonsil tissue lysate

Developed using the ECL technique.

Western blot - Anti-FRRS1 antibody (ab121538)

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

---

### **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors