

Product datasheet

Anti-Fascin-3 antibody ab91484

2 图像

概述

|       |                                                                                                        |
|-------|--------------------------------------------------------------------------------------------------------|
| 产品名称  | Anti-Fascin-3抗体                                                                                        |
| 描述    | 兔多克隆抗体to Fascin-3                                                                                      |
| 经测试应用 | 适用于: WB, ELISA                                                                                         |
| 种属反应性 | 与反应: Mouse, Human                                                                                      |
| 免疫原   | Synthetic peptide selected from the N-terminal region of Human Fascin-3 conjugated to KLH (NP_065102). |
| 阳性对照  | NCI-H460 cell line lysate and mouse testis tissue lysate.                                              |

性能

|      |                                                                                                         |
|------|---------------------------------------------------------------------------------------------------------|
| 形式   | Liquid                                                                                                  |
| 存放说明 | Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.                   |
| 存储溶液 | Preservative: 0.09% Sodium Azide<br>Constituents: PBS                                                   |
| 纯度   | Ammonium Sulphate Precipitation                                                                         |
| 纯化说明 | ab91484 is prepared by Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS. |
| 克隆   | 多克隆                                                                                                     |
| 同种型  | IgG                                                                                                     |

应用

Our [Abpromise guarantee](#) covers the use of **ab91484** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

| 应用    | Ab评论           | 说明 |
|-------|----------------|----|
| WB    |                |    |
| ELISA |                |    |
| 应用说明  | ELISA: 1/1000. |    |

WB: 1/50 - 1/100. Predicted molecular weight: 57 kDa.


Not yet tested in other applications.

Optimal dilutions/concentrations should be determined by the end user.

## 靶标


|       |                                    |
|-------|------------------------------------|
| 功能    | Acts as an actin bundling protein. |
| 组织特异性 | Expressed in testis.               |
| 序列相似性 | Belongs to the fascin family.      |
| 细胞定位  | Cytoplasm > cytoskeleton.          |

## 图片



Anti-Fascin-3 antibody (ab91484) at 1/50 dilution + NCI-H460 cell line lysate at 35  $\mu$ g

**Predicted band size : 57 kDa**



Anti-Fascin-3 antibody (ab91484) at 1/50 dilution + Mouse testis tissue lysate at 35  $\mu$ g

**Predicted band size : 57 kDa**

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish

- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors